

AQUACLEAN mini

MAKE YOUR SKIN YOUNG

Pore
Cleansing

Facial
Lifting

Solution
Delivery

The combination of
3 functions for
complete
skin care

With valve
finely adjustable
amount



AQUACLEAN mini

MAKE YOUR SKIN YOUNG

3-in-1

AquaClean mini consist of 3-in-1 system that combines three different technologies including the aqua clean system using **vacuum**, **electroporation** and **double lifting**. 3-in-1 system provides all facial care including cleansing, solution penetration and lifting to a single device.



① Aqua Cleansing

Valve type
Vacuum & Flow rate control



② Multi-polarity
Electroporation



③ Double Lifting

Cleansing & Exfoliation
by vacuum



Removal effect of sebum,
white head, black head



Effective skin cleansing



Various solutions
can be applied



The effects of specialized
facial lifting treatment



Penetration of solution
to the deep dermis



Premium Solution



1 Remove sebum
and oxidized waste
deeper in the pores
by addition of
Lactic Acid.



2 Dissolve oxidized
dead skin and calms
acne, rash, and other
skin problems by
addition of
Salicylic Acid.



3 Supplement
nutrients and deep
moisturizing on
the cleaner skin by
addition of Portulaca
oleracea extract.

AQUACLEAN mini 3-in-1 System

Optimum Multi-functional skin treatment

Hypoallergenic
peeling function



+

Solution
Penetration

+

Double lifting

1 Aqua Clean System



- Hypoallergenic peeling function that removes sebum, dead skin cell and wastes using vacuum.
- Simple treatment procedure without any side-effects.
- Smooth peeling & treatment without stimulation.



Small



Large



Cap

2 Electroporation

- Without damage to the cellular skin structure, micro-current creates a temporary opening on the skin for penetration of nourishing ampoules into deep layers of the skin.
- Adjustable Treatment Level (Level 1 to 10)

3 Double Lifting

- By supplying micro-current to flow through superficial layer of epidermis, lifting effect from stimulated facial muscle is achieved.
- Adjustable Treatment Level (Level 1 to 10)

4 Diamond Peel

- Diamond peel head tip for improving skin elasticity of the cheeks, neck, neckline, etc.
- Improving skin complexion.



Head Tip

④ Diamond Peel



AquaClean Mini has the option to add another facial care that is suitable for medical skin clinics and hospital centers, making it a **4-in-1 system** that combines aqua clean system using vacuum, electroporation, double lifting and diamond peel.

SPECIFICATION

Applicators	Aqua cleansing / Electroporation Double Lifting
Input Voltage	AC 100-120V/220-240V, 50/60Hz
Weight(kg)	17kg
Dimension(WxDxH mm)	430 x 380 x 380
Exclusively distributed by	Spec Brite Group (GCC and Middle East)



Lebanon
+961 9 215709
+961 71 265 818

Abu Dhabi
+971 2 6275926
+971 2 6267370
P.O. Box 48427

Dubai
+971 4 3884549
+971 4 3881152
P.O. Box 114467

info@specbritegroup.com | www.specbritegroup.com